

# CLIENT BULLETIN



*December 1, 2004*

*Donald R. Simpson, M.D, PhD.*

Dear Valued Customer:

Pactox will have a change in blood transport containers. Pactox will be discontinuing the foam tube holder insert gradually over the next 60 days. The foam tube holder insert will be replaced with the Aqua-Pak Blood Kit. The Aqua-Pak Blood Kit is made of absorbent material and has 4 and/or 6 tube holder pockets. This change is due to an OSHA regulation. Instructions for handling the Aqua-Pak Blood Kit is as follows:

1. Insert one tube per pocket
2. Place Aqua-Pak bundle in a leak-proof **Ziplock** bag
3. Place the finished Aqua-Pak bundle into the provided Pactox lab box

If you have further questions or concerns, please do not hesitate to contact us. We can be reached at our toll free number: 800-328-6942. Ask for Jose Guzman at Ext. 413. Email communication is also an option. Jose Guzman email address is [jguzman@pactox.com](mailto:jguzman@pactox.com).

---

Donald R. Simpson, MD, PhD  
Medical Director  
Pacific Toxicology Laboratories

---

Kingsley R. Labrosse, PhD., DABCC(T)  
Laboratory Director  
Pacific Toxicology Laboratories