

CLIENT BULLETIN

July 2018

Amphetamine/ Methamphetamine

Dear Valued Client,

Effective 7/2/2018, the quantitative cut-off level for Urine Amphetamine/Methamphetamine has been lowered to 100 ng/mL. The purpose of this change is to better assist the physician to determine patient compliance.

Test Code	Test Name	Method	Quantitation Level	Turnaround Time
602	Amphetamine/ Methamphetamine, Urine	LCMS/MS	100 ng/mL	1-2 Days

Our technical staff is available for consultation at your convenience. You may also contact your sales representative for additional details.



Donald Simpson, MD
Medical Director
Pacific Toxicology Laboratories